

Automated apolipoprotein quantification in blood

using ApoEdge and targeted mass spectrometry

A robust, reproducible and scalable workflow for absolute proteomic quantification.

Aleksandra Pel¹, Rafał Szewczyk¹ and David Kotoł²

1. LabExperts sp z o.o., Gdańsk, Poland
2. ProteomEdge AB, Stockholm, Sweden

Abstract

This study evaluates an automated sample preparation workflow for the ApoEdge™ (ProteomEdge™ AB, Stockholm, Sweden) targeted LC–MS/MS multiplex panel, designed for the absolute quantification of human apolipoproteins in serum. Serum samples were prepared, including digestion and solid-phase extraction, using a SOLO™ Automated Liquid Handler (Hudson Lab Automation, Springfield, NJ, USA), and subsequently analyzed on a QTRAP® 6500+ mass spectrometer (SCIEX, Framingham, MA, USA). Method repeatability, assessed across five independent replicates, demonstrated high precision, with coefficients of variation (CV) consistently remaining below 10%. By eliminating manual pipetting errors and significantly reducing hands-on time, this robust automated platform provides a standardized, scalable solution ideally suited for high-throughput clinical proteomics and large-scale cohort studies.

Introduction

Apolipoproteins are key components of lipid metabolism and important biomarkers in cardiovascular diseases, metabolic disorders, and other clinical conditions. Reliable and reproducible quantification of apolipoproteins is therefore essential in both clinical and proteomics research.

ApoEdge (ProteomEdge AB, Stockholm, Sweden), a panel of heavy-labeled internal protein standards, is a targeted LC–MS/MS tool developed for multiplex absolute quantification of human apolipoproteins in plasma or serum. The panel contains heavy-labeled Quantitative Recombinant Protein Standards (qRePS™) for all human apolipoproteins, enabling accurate multiplex absolute quantification. All ProteomEdge products are delivered pre-aliquoted and dried




Figure 1: qRePS are recombinantly expressed protein fragments of around 100 amino acids that undergo simultaneous processing as their corresponding endogenous proteins. Samples are spiked on top of the qRePS internal standards as the first step of sample preparation providing a strategy to control analytical variability and enable multiplex absolute protein quantification.

in a standardized 96-well plate format with individual spike-ins adjusted for digestion together with one microliter of plasma or serum and can be easily implemented in automated workflows (Figure 1).

Automation of sample preparation has become increasingly important in LC–MS/MS workflows, particularly in high-throughput proteomics and clinical research applications. Automated liquid handling minimizes manual intervention, reduces variability associated with pipetting and sample processing, improves reproducibility, and enhances workflow standardization across large sample sets. In addition, automation supports efficient processing in 96-well plate formats, enabling scalable and robust sample preparation for quantitative proteomic analyses. In this work, serum sample preparation for the ApoEdge workflow was performed using the SOLO Automated Liquid Handler (Hudson Lab Automation, Springfield, NJ, USA), allowing automated processing directly in qRePS-coated 96-well plates prior to LC–MS/MS analysis using a QTRAP 6500+ mass spectrometer (SCIEX, Framingham, MA, USA).

Materials and Methods

Reagents and Standards

For the automated digestion procedure, the following reagents were used: 0.3% RapiGest prepared in ultrapure water, 30 mM DTT (Dithiothreitol) prepared in 1x PBS (Phosphate buffered saline), 200 mM CAA (2-Chloroacetamide) prepared in 1x PBS and trypsin (0.1 µg/µL) prepared in 1x PBS. For the solid-phase extraction (SPE) procedure, MeOH (Methanol), Solution A (0.5% acetic acid in water), and Solution B (0.5% acetic acid, 80% acetonitrile, 20% water) were used. All solvents were of LC–MS grade.

Human serum was obtained from UTAK Laboratories, Inc. (Valencia, CA, USA) and used as the biological matrix. The ApoEdge panel of qRePS targeting all human apolipoproteins was used for targeted LC–MS/MS absolute quantification.

Sample Preparation

The development of the automated sample preparation workflow was based on the ProteomEdge protocol “Digestion RapiGest” (https://proteomedge.com/wp-content/uploads/2025/12/Digestion_RapiGest_Rev4_Dec2025.pdf).

Automated sample preparation was performed using the SOLO Automated Liquid Handler according to the optimized layout (Figure 2) and sample preparation workflow (Figure 3). First, serum samples were transferred into the ApoEdge plate and subsequently subjected to denaturation, reduction, alkylation, and enzymatic digestion. Subsequently, layout of reagents and consumables were set-up according to the layout (Figure 4) and SPE was performed on the digested samples according to the protocol (Figure 5). After SPE, the samples were concentrated by evaporation using a Concentrator Plus (Eppendorf, Hamburg, Germany), operated in V-HV (vacuum – high vapor) mode at 45°C for 2 hours. The sample preparation workflow is presented in the schemes below.

The performance of both the automated digestion and SPE procedures was evaluated in terms of repeatability. For this purpose, five independent replicates were prepared and processed under identical conditions.



Figure 2: Layout of plates, tips, and heating block on the Hudson SOLO Automated Liquid Handler during the digestion procedure.




Figure 3: Overview of the automated sample preparation workflow for the Hudson SOLO Automated Liquid Handler with all the steps required during the digestion procedure.



Figure 4: Layout of plates, tips, and heating block on the Hudson SOLO Automated Liquid Handler during the solid-phase extraction procedure.




Figure 5: Overview of the automated sample preparation workflow for the Hudson SOLO Automated Liquid Handler with all the steps required during the solid-phase extraction procedure.

Chromatographic and Mass Spectrometric Conditions

Chromatographic separation was performed in reversed-phase mode using a Nexera LC-40 System (Shimadzu Corporation, Kyoto, Japan) equipped with an ACQUITY UPLC® CSH C18 Column (Waters Corporation, Milford, MA, USA) (130 Å, 1.7 µm, 2.1 × 100 mm) and an ACQUITY UPLC CSH C18 VanGuard™ Pre-Column (Waters Corporation, Milford, MA, USA) (130 Å, 1.7 µm, 2.1 × 5 mm). A gradient elution at 500µL/min was applied using mobile phase A consisting of 0.1% formic acid in water and mobile phase B consisting of 0.1% formic acid in acetonitrile.

Mass spectrometric analysis was performed using the SCIEX QTRAP 6500+ MS/MS System operating in Multiple Reaction Monitoring (MRM) mode with electrospray ionization (ESI). Detailed chromatographic conditions, gradient parameters, mass spectrometer settings, and all MRM transitions are provided in the [Supporting Materials](#).

Results

Data processing method was developed using SCIEX OS Software with the Analytics Module and MQ4 Integration Algorithm (SCIEX, Framingham, MA, USA). To assess repeatability, the coefficient of variation (CV) was calculated based on quantitative results for each peptide. Due to the large number of results, one representative peptide per protein was presented. The CV across five independent replicates did not exceed 7% with median CV of 2.9% (Table 1), confirming high method repeatability indicating that the automated sample preparation, MS analysis and absolute quantification can be reliably applied under laboratory conditions and across wide dynamic range of serum protein concentrations (Figure 6).

For Apolipoprotein F and Apolipoprotein A-V, only signals corresponding to heavy-labeled peptides were detected. Therefore, quantitative results could not be calculated, and these proteins were not included in the table. In addition, Apolipoprotein A-V is still under method optimization.

Table 1: Absolute serum protein concentrations for each biological replicate and their corresponding coefficients of variation (CV).

Protein_Peptide	Serum Concentration [pmol/µL]					Average	SD	CV [%]
	Sample 1	Sample 2	Sample 3	Sample 4	Sample 5			
APOA2_SPELQAEAK	55.067	57.332	56.042	55.187	54.432	55.612	1.120	2.0
APOA1_VEPLRAELQEGAR	31.225	31.288	31.199	31.926	30.402	31.208	0.541	1.7
APOC1_EFGNTLEDK	11.677	13.829	13.768	13.865	13.247	13.277	0.929	7.0
APOC3_DALSSVQESQVAQQAR	5.934	6.115	5.848	5.968	5.588	5.891	0.195	3.3
APOB_ITENDIQIALDDAK	5.522	6.037	5.628	5.942	5.755	5.777	0.214	3.7
CLU_ASSIIDELFQDR	2.261	2.299	2.256	2.246	2.255	2.263	0.021	0.9
APOH_VCPFAGILENGAVR	1.989	2.068	1.987	2.039	1.986	2.014	0.038	1.9
APOA4_SLAELGGHLDQQVEEFR	1.298	1.317	1.314	1.273	1.320	1.304	0.020	1.5
APOM_AFLLTTPR	0.941	0.980	0.935	0.952	0.929	0.947	0.020	2.1
APOE_EQVAEVR	0.753	0.775	0.731	0.731	0.722	0.743	0.022	2.9
APOL1_VTEPISAESGEQVER	0.695	0.693	0.690	0.680	0.664	0.685	0.013	1.9
APOD_VLNQELR	0.362	0.394	0.383	0.392	0.371	0.380	0.014	3.6
APOC4_AWFLESK	0.142	0.152	0.156	0.153	0.141	0.149	0.007	4.5
LPA_GTYSTTVTGR	0.124	0.112	0.115	0.123	0.115	0.118	0.006	4.7
APOC2_STAAMSTYTGIFTDQVLSVLK	0.066	0.065	0.063	0.058	0.068	0.064	0.004	5.8




Figure 6: Absolute serum concentrations of the fifteen different apolipoproteins quantified in 5 biological replicates with the CVs distribution and median CV of 2.9%.

Good chromatographic separation and highly stable retention times were achieved across all five biological replicates, demonstrating excellent reproducibility of the LC-MS/MS method under analytical flow conditions (Figure 7). Overlay of the total ion chromatograms (TICs) showed consistent peptide elution profiles with minimal retention time drift between runs, indicating robust chromatographic performance and stable instrument operation throughout the analysis. The combination of automated sample preparation, reproducible chromatographic separation, and targeted MRM detection enabled reliable multiplex quantification of apolipoproteins over a broad dynamic concentration range. The robustness and reproducibility of the workflow support its suitability for routine high-throughput applications, enabling processing and analysis of up to 96 samples per day in clinical proteomics and large-scale cohort studies.




Figure 7: Overlay of total ion chromatograms (TIC) for the five independent replicates.

Conclusion

- The automated sample preparation in combination with robust MS readout and implementation of internal standards early in the workflow showed high repeatability, as demonstrated by five independent replicates with a maximum CV not exceeding 10%, confirming that the method is robust and reliable for quantitative applications.
- The use of the Hudson SOLO Automated Liquid Handler offers significant advantages over manual sample processing, including reduced pipetting and dilution errors, minimized operator-dependent variability, and improved consistency both within and between plates. In addition, automation substantially decreases hands-on time improving laboratory efficiency and enabling more effective allocation of resources.
- The automated workflow together with ready-to-use panels supports high-throughput sample processing and is therefore well suited for large-scale cohort studies, enabling simultaneous preparation of multiple samples with high consistency.
- Standardization and integrated quality control throughout the complete workflow enhance method transferability, reproducibility, and scalability, which are critical for clinical and proteomics applications. Reproducible workflows implementing panels of internal protein standards, such as ApoEdge, are key enablers of reproducible and harmonized sample preparation, ensuring robust and reliable quantitative results.